# EDIT - bug #7002

## Delete of contig sequence not possible

10/04/2017 10:59 AM - Patrick Plitzner

Status: New Start date:

Priority: New Due date:

Assignee: Katja Luther % Done: 0%

Category: cdmlib Estimated time: 0:00 hour

Target version: Unassigned CDM tickets

Severity: critical Found in Version:

### Description

Deletion in specimen hierarchy editor is not possible because of the following message.

An object of eu.etaxonomy.cdm.model.molecular.DnaSample with ID 9381 is referencing the object

It is also not possible to configure the deletion process so that it works.

#### History

#### #1 - 10/04/2017 11:39 AM - Katja Luther

we decided that it is not possible to delete media from media specimen and from other places where they are directly associated to an cdm object.

you have to delete it from the dna sample then you can delete it in bulk editor

### #2 - 10/04/2017 11:54 AM - Patrick Plitzner

Katja Luther wrote:

we decided that it is not possible to delete media from media specimen and from other places where they are directly associated to an cdm object.

you have to delete it from the dna sample then you can delete it in bulk editor

No, sorry, I meant the CDM entity Sequence with "contig sequence" I just build a simple hierarchy FU->specimen->DnaSample->Contig Sequence. No special info or details added. And then I try to delete the contig sequence in the specimen editor via context menu.

03/20/2024 1/1